

Tutorial #X: Preliminary Report

Barry Linkletter

This is the final report produced by your unfortunate predecessor. It describes the enzyme assays for wild-type β -galactosidase with PNP- β -D-Gal. Please use this report to develop your own method for analyzing the remaining data that still needs to be analyzed. Then you may have cake.

This document was produced using the \LaTeX typesetting language with the Tufte-handout document class. Images of proteins were created using UCSF Chimera. Chemical diagrams were made with ChemDoodle and further edited with Affinity Designer.

Part 1: Introduction

This report describes the data analysis for the enzyme assay used to evaluate mutated β -galactosidase enzymes. First we must determine an enzyme concentration that will give reaction rates that are fast enough to be easily followed but not too fast where we cannot collect enough data to determine an initial rate.

The Plate Plan

We will evaluate three different concentrations of enzyme. The plate will be set up with eight rows of differing substrate concentration in 0.100 M phosphate buffer. Three rows will contain no enzyme, three will have an enzyme concentration of 1 nM and there will be two sets of three columns with dilutions of that concentration. In the first trial we will use 2-fold dilutions. The plate plan is sketched in figure 1.

		[Enzyme] /nM											
		0			1.0			0.5			0.25		
[PNP- β -Gal] /mM		1	2	3	4	5	6	7	8	9	10	11	12
0.01	A	○	○	○	○	○	○	○	○	○	○	○	○
0.02	B	○	○	○	○	○	○	○	○	○	○	○	○
0.03	C	○	○	○	○	○	○	○	○	○	○	○	○
0.04	D	○	○	○	○	○	○	○	○	○	○	○	○
0.06	E	○	○	○	○	○	○	○	○	○	○	○	○
0.10	F	○	○	○	○	○	○	○	○	○	○	○	○
0.20	G	○	○	○	○	○	○	○	○	○	○	○	○
0.50	H	○	○	○	○	○	○	○	○	○	○	○	○

Figure 1: Plate plan for determining optimal enzyme concentration for enzyme assays of β -galactosidase with PNP- β -D-Gal

Two 96-well plates were set up. One contained the enzyme in phosphate buffer and the other PNP- β -D-Gal in a 8 % mixture of isopropanol and phosphate buffer. The phosphate buffer had a concentration of 0.100 M and was at pH 7.0. Volumes in each plate were equal. At time zero the contents of substrate plate was transferred to the enzyme plate by a multipipettor. The series of samples were prepared so that the final concentrations in the plate as described in the plate plan were realized.¹ The plate was placed in

¹ Question: After the plates are mixed, what is the amount of isopropanol present in each sample. Do you think it will affect the results? Why did I choose to have a isopropanol as a small fraction of the solution?

a plate reader and readings began 30 seconds after the addition of substrate. Absorbance was measured at 405 nm.²

² Question: What do we know about PNP- β -D-Gal that informed the choice of this wavelength.

Data Collection

Data was exported from the plate reader as separate time/absorbance files names according to column and row. A python script to plot the initial results for each well in a column of the plate was used to produce the contact sheet of data below.




Figure 2: Plots of A_{405} vs time for each well in the plate. Each column represents an experiment at eight different substrate concentrations.

The *Python* code for generating this plot is available via Google Colab.

Observations

Substantial experimental error was observed in the data and the measurements became unreliable when absorbances exceed the value of 2.³

The time frame of the experiment was such that most of the reactions reached completion. For initial rate analysis we must use the absorbance change in the few percent change in substrate concentration as the reaction proceeds. We want the slope to be as close as possible to the initial slope. I wrote a *Python* script that plots data from a single well over a selected time frame, performing a linear fit of the data and calculating the residuals to confirm that we are using a time frame in which the data is mostly linear. We want as much data as possible so would like a longer time frame, but we want to have little or no observable curve in the line as so want a small a time frame as possible after the start of the reaction.

After examining many cells and trying several time frames, I chose a time frame from zero to 3 minutes. This captures about 15 of the 358 data points collected over the hour duration of the experiment.

³ Question: The manufacturer of the plate reader advertised that precise readings could be obtained up to 3 absorbance units. Given what we know about how absorbance is calculated, please explain why the error is much greater at higher absorbance values.




Figure 3: Plots of A_{405} vs time for well 7-B. [Enzyme] = 0.5 nM, [PNP- β -D-Gal] = 0.02 mM. Three time frames are shown: the full 60 min, 12 min, and 4 min. The residuals for the line fit of the 4 min time frame are shown in the lower left plot. No curvature was detectable within the range of the experimental error.

The *Python* code for generating this plot is available via Google Colab.

The short time frame was required because the highest enzyme concentration resulted in reactions so rapid that substantial change in substrate concentration occurred in the first few moments. For wells where the substrate concentration was at or below the K_M value, we saw significant curvature at longer time frames. The diluted enzyme samples were linear over longer time frames and so more data could have been used. I will use the

short time frame discussed above to be consistent but future experiments with less enzyme may be able to use more data points before curvature is observed. Hopefully my unfortunate successor will find this information useful. The cake is a lie!

There is significant hydrolysis of the PNP- β -D-Gal in the absence of enzyme catalyst. It is a very reactive molecule and at pH 7.0 the rate of hydrolysis is significant. It will be necessary to subtract the background rate from catalyzed rates to ensure accuracy.

Data Analysis

Based on the observations above I will create a *Python* script to calculate the initial slope of absorbance change over time in each cell. We can calculate the rate of appearance of product using the molar extinction coefficient for phenolate anion and the pK_a value for the phenol. The curve fit function will also report a standard deviation for the slope based on the random error in the data and the number of points. We will use this value as a measure of the precision of our data. The script will output the rate of reaction and the standard error for all 96 cells and save that data in a csv file. The python code is available as an interactive *Python* notebook via Google Colab.